

Fast, accurate detection helps stop the spread of SARS-CoV-2, Influenza virus, and RSV.

Solutions Powered by Nano-Check Products

Nano-Ditech Corp. NJ, USA

Tel: 1-609-409-3300, Info@nanoditech.com, www.nanoditech.com

Nano-Check™ COVID-19 Antigen Test


CLIA 
WAIVED

FDA
CLEARED

- ✓ FDA 510(k) Cleared (k231187)
- ✓ FDA CLIA Waived (CW240009)
- ✓ Anterior Nasal Swab
- ✓ Results in 15 minutes
- ✓ No Extra Steps
- ✓ Relative Sensitivity: 84.8%
- ✓ Relative Specificity: 99.7%
- ✓ Made in USA



► Easy/Simple To Use



Contents

- 20 Test Devices
- 20 Reagent tubes with extraction buffer (or ampules)
- 20 Sample collection swabs
- 1 Positive control swab
- 1 Negative control swab
- 1 Instructions for Use
- 1 Quick Reference Instruction

► Performance Characteristics

SARS-CoV-2 Variant	LoD
Wild: USA-WA1/2020	7.0×10^2 TCID 50/mL
Omicron: USA/MD-HP20874/2021	1.95×10^2 TCID 50/mL
Delta: USA/PHC658/2021	5.21×10^2 TCID 50/mL
WHO (NIBSC 21/368)	44 IU/swab

- For additional information regarding our diverse range of diagnostic tests, kindly visit our website at www.nanoditech.com
Tel: 1-609-409-3300
Info@nanoditech.com

Refer to IFU for additional LOD / Analytical Reactivity information

Ordering Information
Cat. No.: MD8147CW (20 Test)

ND
NANO-DITECH CORPORATION

Nano-Check™ RSV Test


CLIA 
WAIVED

FDA
CLEARED

- ✓ FDA 510(k) Cleared (k240280)
- ✓ FDA CLIA Waived (CW240005)
- ✓ Anterior Nasal Swab
- ✓ Infant/Child: 6 months - 6 years old
- ✓ Adult Senior: ≥ 60 years old
- ✓ Results in 15 minutes
- ✓ No Extra Steps
- ✓ Relative Sensitivity: 83.8%
- ✓ Relative Specificity: 99.6%
- ✓ Made in USA



► Easy/Simple To Use



Contents

- 25 Test Devices
- 25 Ampules
- 25 Sample collection swabs
- 1 Positive control swab
- 1 Negative control swab
- 1 Instructions for Use
- 1 Quick Reference Instruction

► Performance Characteristics

RSV subtype	Strain	LoD
A	Long	1.6×10^3 TCID 50/mL
	A-2	2.0×10^4 TCID 50/mL
B	18537	4.0×10^3 TCID 50/mL
	B1	3.5×10^3 TCID 50/mL

- For additional information regarding our diverse range of diagnostic tests, kindly visit our website at www.nanoditech.com
Tel: 1-609-409-3300
Info@nanoditech.com

Refer to IFU for additional LOD / Analytical Reactivity information

Ordering Information
Cat. No.: MD8152 (25 Test)


ND
NANO-DITECH CORPORATION

Nano-Check™ Influenza+COVID-19 Dual Test

- ✓ FDA 510(k) Cleared (K243561)
- ✓ FDA CLIA Waived (CW240029)
- ✓ Anterior Nasal Swab
- ✓ Evaluation conducted with an adequate number of positive samples
- ✓ Results in 15 minutes
- ✓ No Extra Steps
- ✓ Made in USA



► Easy/Simple To Use



Contents

- 25 Test Devices
- 25 Ampules
- 25 Sample collection swabs
- 1 Positive control swab
- 1 Negative control swab
- 1 Instructions for Use
- 1 Quick Reference Instruction

► Performance Characteristics

SARS-CoV-2	Comparator RT-PCR		Total
	Positive	Negative	
Positive	232	3	235
Negative	33	1701	1734
Total	265	1704	1969

Positive Percent Agreement = 87.6%
 Negative Percent Agreement = 99.8%

Flu A	Comparator RT-PCR		Total
	Positive	Negative	
Positive	417	6	423
Negative	63	1483	1546
Total	480	1489	1969

Positive Percent Agreement = 86.9%
 Negative Percent Agreement = 99.6%

Flu B	Comparator RT-PCR		Total
	Positive	Negative	
Positive	99	5	104
Negative	15	1850	1865
Total	114	1855	1969

Positive Percent Agreement = 86.8%
 Negative Percent Agreement = 99.7%

Refer to IFU for additional LOD/Analytical Reactivity/Clinical Performance Information

Ordering Information
Cat. No.: MD8154 (25 Test)


► For additional information regarding our diverse range of diagnostic tests, kindly visit our website at www.nanoditech.com Tel: 1-609-409-3300, Info@nanoditech.com

Nano-Check™ Influenza A+B Test


CLIA 
WAIVED

FDA
CLEARED

- ✓ FDA 510(k) Cleared (K252283)
- ✓ FDA CLIA Waived (CW250012)
- ✓ Anterior Nasal Swab
- ✓ Evaluation conducted with an adequate number of positive samples
- ✓ Results in 15 minutes
- ✓ No Extra Steps
- ✓ Made in USA



▶ Easy/Simple To Use



Contents

- 25 Test Devices
- 25 Ampules
- 25 Sample collection swabs
- 1 Positive control swab
- 1 Negative control swab
- 1 Instructions for Use
- 1 Quick Reference Instruction

▶ Performance Characteristics

Flu A	Comparator RT-PCR		Total
	Positive	Negative	
Positive	417	6	423
Negative	63	1483	1546
Total	480	1489	1969

Positive Percent Agreement = 86.9%
Negative Percent Agreement = 99.6%

Flu B	Comparator RT-PCR		Total
	Positive	Negative	
Positive	99	5	104
Negative	15	1850	1865
Total	114	1855	1969

Positive Percent Agreement = 86.8%
Negative Percent Agreement = 99.7%

Refer to IFU for additional LOD/Analytical Reactivity/Clinical Performance Information

Ordering Information
Cat. No.: MD8153(25 Test)

- ▶ For additional information regarding our diverse range of diagnostic tests, kindly visit our website at www.nanoditech.com Tel: 1-609-409-3300, Info@nanoditech.com